

# Omnipod® System and Dexcom or Eversense CGM Data Integrated in Insulet Provided Glooko!



**IMPORTANT:** Before making any changes to a user's PDM settings that may affect insulin delivery, please review current settings and history summaries within the user's PDM.



**How?** Patient must connect their Dexcom or Eversense account to their Glooko account. Once connected, their Omnipod® System and CGM data will be available in Glooko.

**Patient can connect by following these 3 simple steps:**

- 1** Log into Glooko account on [my.glooko.com](https://my.glooko.com)
- 2** Click **Settings** in the top right
- 3** Scroll to the bottom and click **Connect Dexcom Account** or **Connect Eversense Account**



**Don't forget to enter your clinic's ProConnect code into your account!** That way, you can upload your Omnipod® System at home and your healthcare provider will be able to see your latest data.

**If you have any questions about Glooko, please contact Glooko Support:**

**Call:** 1-800-206-6601 **Text:** 650-720-5310 **Email:** [support@glooko.com](mailto:support@glooko.com)

---

**For more information on Omnipod®, contact Insulet Customer Support:**

**Call:** 1-800-591-3455 **Visit:** [omnipod.com](http://omnipod.com)



**Insulet Corporation** 100 Nagog Park, Acton, MA 01720 • 800-591-3455

© 2021 Insulet Corporation. Omnipod, the Omnipod logo, and Simplify Life are trademarks or registered trademarks of Insulet Corporation. All rights reserved. Glooko is a trademark of Glooko, Inc. and used with permission. Dexcom is a registered trademark of Dexcom, Inc. and used with permission. All other trademarks are the property of their respective owners. The use of third party trademarks does not constitute an endorsement or imply a relationship or other affiliation. INS-ODS-03-2021-00126 V1.0