

A woman with dark hair, smiling, wearing a yellow sleeveless top and yellow shorts. She has a white insulin pump on her waist. She is posing with one leg raised. The background is a solid yellow color with a large, rounded yellow shape in the center.

MAKE DIABETES A SMALLER PART OF LIFE

For people with type 1 diabetes,
ages 2 years and older.

omnipod[®]
5
automated insulin
delivery system

LIFE GETS BIGGER.



Pod shown without
necessary adhesive.

when diabetes management becomes smaller.

Multiple daily injection-free. It's automated insulin delivery, simplified.

Dealing with unexpected highs and lows? **Not fine.**

Waking up in the middle of the night? **Not necessary.**

Interrupting the moment for an injection? **No more.**

Status quo is overrated. Good enough is not good enough. You deserve more.

Today is the day to try something new.

GET TANGLED UP IN LIFE, NOT TUBES.

It's time for diabetes to take a step back.

Omnipod® 5 lets you work, exercise, shower, and swim on your terms, without a second thought. Omnipod 5 is also a simple choice for kids—allowing them to worry less about diabetes and focus on what matters most in their exciting lives. They can explore, learn, climb, run, play, and splash without a second thought. Plus, this sleek, waterproof*, and tube-free wearable Pod is discreet and easy to conceal under clothing—so diabetes is only part of the conversation if you want it to be.

Omnipod 5 automatically delivers insulin to help protect against highs and lows, day and night, so you can spend more time living life and less time thinking about diabetes.^{1,2}

*The Pod has a waterproof IP28 rating for up to 7.6 metres (25 feet) for 60 minutes. The Controller is not waterproof.



**THERE
ARE LOTS
OF PUMPS.
THERE'S
ONLY ONE
OMNIPOD 5.**




**THINK
ABOUT A LIFE
BEYOND
MDI!**

**Omnipod 5 users agree
that it is easy to use***

Pod shown without necessary adhesive.
*Insulet Data on File. OP5-003 Clinical Study Report. 2024. A 13-week randomized controlled trial conducted among 194 adults (age 18-70 years) with type 1 diabetes in France and the U.S., comparing the safety and effectiveness of the Omnipod 5 Automated Insulin Delivery (AID) System versus conventional non-AID pump therapy and CGM (control). At 13-weeks of Omnipod 5 use, 130 participants were asked via a System Usability Survey (SUS) if the system was easy to use. 71.5% of participants strongly agreed and 19.2% agreed.

3

SIMPLE PARTS.



Sensor*

Pod shown without necessary adhesive.

Omnipod 5 Controller

Take control with the Omnipod 5 handheld device.


Pod

Tubeless, wearable, and waterproof* with built-in SmartAdjust™ technology—a predictive insulin algorithm—the Pod sits right on your body and helps to keep glucose in range.^{1,2}

Your Sensor

In constant communication, your sensor* sends glucose values and trends to the Pod to help the System make automated insulin delivery decisions. Please visit omnipod.com for a list of compatible sensors.

Learn more about Omnipod 5 and have any questions answered live at a Pod 101 webinar.



Scan to register for a session at omnipod.com.

*Sensor prescribed and sold separately. Compatible sensor options may vary by country.

¹The Pod has a waterproof IP28 rating for up to 7.6 metres (25 feet) for 60 minutes. The Controller is not waterproof.

**3 EASY STEPS,
EVERY 3 DAYS***

1.



2.



3.



1. Fill

Fill the Pod with up to 200 units of rapid-acting insulin and pair with Controller.

2. Place

Place the Pod almost anywhere you would give an insulin injection. The cannula (a small tube) inserts automatically with the push of a button.

3. Tap

Use your Controller to control your insulin delivery with a few taps on a touch screen.

ONE POD REPLACES UP TO 14[†] INJECTIONS

Pod shown without necessary adhesive.


*See User Guide for full instructions.

[†]14 injections/3 days based on people with T1D on MDI taking ≥ 3 bolus and 1-2 basal injections/day multiplied by 3 days. Chiang JL, et al. Type 1 Diabetes Through the Life Span: A Position Statement of the American Diabetes Association. *Diabetes Care*. 2014;37(7):2034-2054.

ADJUSTING— SO YOU DON'T HAVE TO.

SmartAdjust™ technology is the brains behind the Omnipod 5 System, helping to protect against highs and lows.^{1,2} The System predicts where glucose values will be in 60 minutes and automatically increases, decreases, or pauses insulin every 5 minutes based on your sensor value and trend. Bolus doses can be given for meals using the Omnipod 5 Controller.

SmartAdjust™ technology: How it works



The Pod has a waterproof IP28 rating for up to 7.6 metres (25 feet) for 60 minutes. The Controller is not waterproof.



SAY HELLO TO RESTFUL NIGHTS.

Omnipod 5 automatically adjusts your insulin levels every 5 minutes, even overnight, to help keep glucose in range.^{1,2} It's the system that never sleeps, so you can start dreaming big.

MEALTIME MATH, DONE FOR YOU.

The SmartBolus Calculator makes it simple to determine the right amount of insulin for meals and snacks. It uses your sensor value and trend along with carbohydrate amount to calculate a precise mealtime dose—no need for guessing.



TECHNOLOGY, MADE SIMPLE FOR YOU.



GET GOING WITH THE ACTIVITY FEATURE.

When turned on, the Activity Feature reduces insulin delivery for times when glucose typically goes low, like exercising. Whether hitting the gym, enjoying a yoga class, or participating in your favourite sport, Omnipod 5 makes the balancing act easier and protecting against lows simple.



CONTROLLED BY A HANDHELD DEVICE.

Take full control of the System at any time with just a few taps. With the Omnipod 5 Controller, Bolusing for a meal, enabling the Activity Feature, and even changing the Pod is so simple.

OMNIPOD 5
IS SMALL.
YET ITS
IMPACT IS

HUGE

Omnipod 5 users around the world trust Omnipod 5 because it helps take their minds off insulin management, giving them more control over both their lives and their diabetes.



“It’s good when he wakes up in range. It helps him feel good and start his day on the right foot.”

LIZ D.
Mom to Sponsored Podder®, Christian



“Omnipod 5 working its magic in the background means she has more time to have fun and worry less about diabetes.”

REBECCA R.
Mom to Podder®, Ava R.



“Some of the decisions I have to make daily regarding diabetes, the device handles it all in the background.”

ALVIN M.
Omnipod 5 user



“The Omnipod 5 has dramatically taken the focus from my diabetes to me focusing on my life, my job, my friends, and my relationships.”

MORGAN M.
Omnipod 5 user

TODAY'S THE DAY TO GAIN PEACE OF MIND.

Ready to experience life with fewer diabetes interruptions? Kicking multiple daily injections and bulky pumps out of your life means more freedom, greater peace of mind, and more time for you and your loved ones.

Order a Free Demo Pod* Today!

Curious to know how it feels to wear one of our tubeless, waterproof† Pods?

Order a free Pod Experience Kit* today and you'll receive a real-feel, real-size, real-weight demo Pod*, without the insulin.

Got questions? We're happy to help!

Our Omnipod Specialists are here to support you on your journey to Pod Therapy.

Call us at **1-855-POD-INFO (1-855-763-4636)**
or email omnipodcanada@insulet.ca.

ORDER A FREE DEMO POD TODAY AT OMNIPOD.COM



*The demo Pod is a needle-free Pod that does not deliver insulin. Controller not included.

†The Pod has a waterproof IP28 rating for up to 7.6 metres (25 feet) for 60 minutes. The Controller is not waterproof.

*Your call may be recorded for quality monitoring and training purposes.

1. Brown S, et al. *Diabetes Care*. 2021;44:1630-1640. Prospective pivotal trial in 240 participants with T1D aged 6 - 70 yrs. Study included a 14-day standard therapy (ST) phase followed by a 3-month Omnipod 5 hybrid closed-loop (HCL) phase. Mean time >10.0 mmol/L or >180mg/dL (12AM-<6AM) as measured by CGM in adults/adolescents and children ST vs. 3-mo Omnipod 5: 32.1% vs. 20.7%; 42.2% vs 20.7%, $P<0.0001$, respectively. Mean time >10.0 mmol/L or >180mg/dL (6AM-<12AM) as measured by CGM in adults/adolescents and children ST vs. 3-mo Omnipod 5: 32.6% vs. 26.1%; 46.4% vs 33.4%, $P<0.0001$, respectively. Mean time <3.9 mmol/L or <70 mg/dL (12AM-<6AM) as measured by CGM in adults/adolescents and children ST vs. 3-mo Omnipod 5: 3.64% vs. 1.17%, $P<0.0001$; 2.51% vs. 1.78, $P=0.0456$, respectively. Mean time <3.9 mmol/L or <70 mg/dL (6AM-<12AM) as measured by CGM in adults/adolescents and children ST vs. 3-mo Omnipod 5: 2.64% vs. 1.37%, $P<0.0001$; 2.13% vs. 1.98%, $P=0.2545$, respectively. Mean time in range (3.9-10.0 mmol/L or 70-180mg/dL) in adults/adolescents as measured by CGM: ST = 64.7%, 3-mo Omnipod 5 = 73.9%, $P<0.0001$. Mean time in range (3.9-10.0 mmol/L or 70-180mg/dL) in children as measured by CGM: ST = 52.5%, 3-mo Omnipod 5 = 68.0%, $P<0.0001$.

2. Sherr J, et al. *Diabetes Care*. 2022; 45:1907-1910. Single-arm multicentre clinical trial in 80 pre-school children (aged 2-5.9 yrs) with T1D. Study included a 14-day standard therapy (ST) phase followed by a 3-month AID phase with Omnipod 5 system. Mean time >10.0 mmol/L or >180mg/dL (12AM-<6AM) as measured by CGM in children ST vs. 3-mo Omnipod 5: 38.4% vs. 16.9%, $P<0.0001$, respectively. Mean time >10.0 mmol/L or >180mg/dL (6AM-<12AM) as measured by CGM in children ST vs. 3-mo Omnipod 5: 39.7% vs. 33.7%, $P<0.0001$, respectively. Mean time <3.9 mmol/L or <70 mg/dL (12AM-<6AM) as measured by CGM in children ST vs. 3-mo Omnipod 5: 3.41% vs. 2.13%, $P=0.0185$. Mean time <3.9 mmol/L or <70 mg/dL (6AM-<12AM) as measured by CGM in children ST vs. 3-mo Omnipod 5: 3.44% vs. 2.57%, $P=0.0799$. Mean time in range (3.9-10.0 mmol/L or 70-180mg/dL) as measured by CGM in children ST vs. 3-mo Omnipod 5: 57.2% vs 68.1%, $P<0.0001$."

Indications for use

The Omnipod 5 Automated Insulin Delivery System is a single hormone insulin delivery system intended to deliver U-100 insulin subcutaneously for the management of type 1 diabetes in persons aged 2 and older requiring insulin.

The Omnipod 5 System is intended to operate as an automated insulin delivery system when used with compatible Continuous Glucose Monitors (CGM). When in automated mode, the Omnipod 5 system is designed to assist people with type 1 diabetes in achieving glycemic targets set by their healthcare providers. It is intended to modulate (increase, decrease, or pause) insulin delivery to operate within predefined threshold values using current and predicted CGM values to maintain blood glucose at variable target glucose levels, thereby reducing glucose variability. This reduction in variability is intended to lead to a reduction in the frequency, severity, and duration of both hyperglycemia and hypoglycemia.

The Omnipod 5 System can also operate in a manual mode that delivers insulin at set or manually adjusted rates.

The Omnipod 5 System is intended for single patient use. The Omnipod 5 System is indicated for use with NovoLog®/NovoRapid®, Humalog®, Trurapi®/Truvelog®/Insulin aspart Sanofi®, Kirsty®, and Admelog/Insulin lispro Sanofi U-100 insulin.

Warning: SmartAdjust™ technology should NOT be used by anyone under the age of 2 years old. SmartAdjust™ technology should also NOT be used in people who require less than 5 units of insulin per day as the safety of the technology has not been evaluated in this population.

The Omnipod 5 System is NOT recommended for people who are unable to monitor glucose as recommended by their healthcare provider, are unable to maintain contact with their healthcare provider, are unable to use the Omnipod 5 System according to instructions, are taking hydroxyurea as it could lead to falsely elevated CGM values and result in over-delivery of insulin that can lead to severe hypoglycemia, and do NOT have adequate hearing and/or vision to allow recognition of all functions of the Omnipod 5 System, including alerts, alarms, and reminders. Device components including the Pod, CGM transmitter, and CGM sensor must be removed before Magnetic Resonance Imaging (MRI), Computed Tomography (CT) scan, or diathermy treatment. In addition, the Controller and smartphone should be placed outside of the procedure room. Exposure to MRI, CT, or diathermy treatment can damage the components.

Warning: DO NOT start to use the Omnipod 5 System or change settings without adequate training and guidance from a healthcare provider. Initiating and adjusting settings incorrectly can result in over-delivery or under-delivery of insulin, which could lead to hypoglycemia or hyperglycemia.



OMNIPOD 5 SIMPLIFY LIFE®

The Pod has a waterproof IP28 rating for up to 7.6 metres (25 feet) for 60 minutes. The Controller is not waterproof.



Get started with Pod Therapy
today at **omnipod.com**



Customer Care: 1-855-POD-INFO (1-855-763-4636)

©2025 Insulet Corporation. Omnipod, the Omnipod logo, Simplify Life and SmartAdjust are trademarks or registered trademarks of Insulet Corporation. All rights reserved. The Bluetooth® word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc. and any use of such marks by Insulet Corporation is under license. All other trademarks are the property of their respective owners. The use of third-party trademarks does not constitute an endorsement or imply a relationship or other affiliation. Insulet Corporation, 1540 Cornwall Rd, Suite 201, Oakville, ON, L6J 7W5, Canada.

INS-OHS-10-2024-00228V2.0